


SHEAR FRAME

Laboratory equipment for shear experiments with consolidated sands and soils



Completion: January 2010
Design & Fabrication: Dani Breitenstein
ETH - Swiss federal institute of technology
Soil and Terrestrial Environmental Physics (STEP)
Universitätstrasse 16
CH-8092 Zürich

Group Website: www.step.ethz.ch

Revision state			
Version	Date	Document name	Administrator
1.0	20.04.2011	ShearFrame manual_1.0	Dani Breitenstein

Contents:

1. DEVICE AND ACCESSORY	3
1.1 INTRODUCTION	3
1.2 DETAILS OF THE SHEAR BOX.....	4
2. SAFETY AND PRECAUTIONS.....	5
3. SPECIFICATION	6

1. Device and Accessory


1.1 Introduction

The process of shearing is the counter movement of two frames located on each other.

The Shear Box is designed to measure acoustic emission. Therefore the maximal displacement is longer and the shear velocity can be adjusted in a broad range. The shear box can work either in a strain controlled mode with constant speed or in a load controlled mode with constant load. Additionally a piston is mounted to the top frame to consolidate the observed soils or granular matter.


The lower frame is longer (9) and mobile in one horizontal direction while the upper frame (10) has no degree of freedom. The upper frame's bearings are force transducers to allow measuring all process forces in horizontal and vertical directions and correcting the calculated values for wall friction. The consolidation as well as the shear movement is measured and logged with displacement sensors.

All measurements are analog and are recorded with a data logger and displayed on a screen with a GUI programmed in LabView. At the bottom of the lower frame is a porous plate which allows setting suction up to 400mbar to the sample material. Suction matrix potential can be measured at two different depths below the shear plane via tensiometer. There are designated bore holes for connecting acoustic sensors inside of the soil.



1.2 Details of the Shear Box


The lower Frame:



The lower frame is longer than the upper frame. Its bottom plate is a sintered stainless steel plate which is porous with an air entry value of approximately 200mbar (max. suction).

The shear length can be reduced to a classical shear experiment with a separating wall so the area of the two frames is equaled (installed on picture).


Consolidation:




The consolidation stress is generated by air pressure and introduced with a pneumatic cylinder.

The air pressure is adjusted with a pneumatic controller. A 3-way valve allows holding the position on a constant pressure, constant consolidation pressure or releasing pressure.

Horizontal force transducer:



Vertical force transducer:



2. Safety and precautions

Shear Box	<p>Be aware the shear force as well as the consolidation force are very high and can seriously injure fingers, hands or other parts in the shear zone. Never touch the apparatus during an experiment even if movements are very slow. The forces exceed the range of small bruises and can definitely cause injury.</p>
-----------	--

3. Specification

Shear Box:	Upper frame:	
	- Area	150 x 150 mm
	- Height	72 mm
	Lower frame:	
	- Area	400 x 150 mm
	- Height	72 mm
	- Max. suction	200 mbar
	Shearing:	
	- Maximal displacement	250 mm
	- Maximal shear stress	250 kN/m ²
	- Minimal shear velocity	0.1 mm/min
	- Maximal shear velocity	200 mm/min
	Consolidation:	
	- Maximal Pressure	125 kN/m ²
	- Maximal Consolidation	25 mm
Actuator:	PERO	Spindle drive Strain controller GS - MOT A Force controller MEC 16 SCU Amplifier: SCU-Module SC-TF Software Perolog 3
Data logger:	National Instruments	NI 9237 4-channel bridge module Ni cdaq chassis NI 9205 32-channel analog input module
Sensors:	Displacement:	
	- Shear Strain:	Resistance, Novotechnik LWG 300LVDT
	- Dilation	LVDT, Inelta Typ ISDT50-K-2420
	Force:	
	- Horizontal:	Beam load cell, Omega LC703-100
	- Vertical:	Beam load cell, Omega LC111-1K
	Consolidation Pressure:	Keller LOE 2